Gebr. Bräsegger GmbH & Co. KG
Trophäenweg Weg 25
32657 Lemgo - Germany
Telefon +49 (0) 5261 701-400
Telefax +49 (0) 5261 701-580
info@kometmedical.de
www.kometmedical.de

© Gebr. Bräsegger GmbH & Co. KG
1977-1997
Reg. 2022-06

- Instructions for use.**
Sterile single-use saw blades for work on bones
- Instrucciones de empleo.**
Hojas de sierra estériles para hueso, uso único
- Gebrauchsanweisung.**
Sägeblätter für Knochen Einmalgebrauch steril
- Consigli d'uso.**
Lame de scie stériles à usage unique pour le traitement de l'os
- Инструкции за употреба.**
Листове за тримане за кости за еднократна употреба стерилни
- Uputa za uporabu.**
Jednokratni sterilni listovi pile za rezanje kosti
- Návod k obsluhu.**
Plátno list na kosti na jednorázové použití, sterilní
- Brugsanvisning.**
Savklinger till engångsbruk steril
- Gebrauchsanweisung.**
Zagblättern tot einmalig gebrauch steril
- Käyttöohjeet.**
Steriloidut ohuilehdet katusatutavud luu saalehdet
- Käyttöohje.**
Luostahan kertäkäyttökäyttö steril
- Objezct upotrebe.**
Amelice osterizovani, upotreba, amputacijom

- Caution, ref. no. 5.4.4:** Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
- Store in a dry place, ref. no. 5.3.4:** Indicates that the medical device has to be protected from moisture.
- Protect from sunlight, ref. no. 5.3.2:** Indicates that the medical devices has to be protected from sources of light.
- Contains hazardous substances, ref. no. 5.4.10:** Indicates a medical device that contains substances that can be carcinogenic, mutagenic, reprotoxic (CMR), or substances with endocrine disrupting properties.

¹ Evidence has shown that even those of the devices that contain more than 0.1 % weight percent of cobalt (CAS: 7440-48-4) do not pose a risk of carcinogenicity, mutagenicity or toxicity for reproduction (CMR classification 1B). Therefore, no further measures are required.

MD	Medical device, ref. no. 5.7.7: Indicates the item is a medical device.
CE	Regulated by European law: CE marking of conformity
CH REP	Non-regulated symbol: Indicates the authorized representative in Switzerland
UK CA	Regulated by UK law: UKCA marking of conformity ("UK Conformity Assessed"), followed by the conformity assessment body's identification number
Rx Only	Non-regulated symbol: Caution: Federal law (USA) restricts this device to sale by or on the order of an physician
QTY	Non-regulated symbol: Quantity
MAT	Non-regulated symbol: Material used
GS1	Barcode: GS1 format
HBC	Barcode: HBC format
UDI	Non-regulated symbol: leading barcode for GUDID

Material: stainless and acid resistant steel

Intended use:
The saw blades for work on bones are combined with a surgical power system and, in some cases, with a saw template. They are intended for work on bones and body tissue of comparable hardness during surgical interventions with a view to treating or alleviating diseases that require the treatment of bones. They are exclusively used on human patients by surgeons and medically trained staff in hospitals or surgical practices in a surgical environment, provided that the potential benefit outweighs the risk. Work on bone is done by the reciprocating or oscillating movements of the saw blades.

Indication:
Osteotomies as part of surgical orthopedic treatments, oral maxillofacial surgery and traumatology

Contra-indications:
- Existing metal allergies and all known contra-indications applicable to surgical interventions.
- Open reactions of tissue or bone necrosis with the use of a motor-powered saw outweighs the potential benefit.
- Motor-powered saws should not be used if there are other contra-indications applicable to orthopedic, surgical interventions.
- Operations which require a certain bone quality or that pose a risk of damaging important tissue structures.

Safety instructions:
- Single-use only. Do not re-use.
- Observe the instructions indicated on the label.
- Before start up, make sure that the instrument is securely locked in place.
- Insert the instrument into the hand piece until it locks.
- The instrument is sterile packed.
- Perform an idle run to ensure the secure fit of the saw blade in the power unit

Warnings:
- Not authorized for work on metal. Do not allow the instrument to come in contact with metal during use.
- Avoid excessive contact pressure.
- Make sure to provide sufficient cooling and irrigation.
- Do not reuse instruments that are damaged in any way.
- During use (in a saw template), prevent jamming or bending of the saw blade. Do not use for leverage.

Residual risks and side effects:
Contamination due to injury with infection to the patient, user or a third person as well as allergic/ toxic reactions of tissue or bone necrosis as well as irreversible damage to the surrounding tissue can be caused by excessive heat generation, soiled/contaminated, non-biocompatible or broken/ damaged instruments, by the absence of indication of the instruments or damage to the tissue. Instruments, pointed edges or sharp blades or by the instruments coming in contact with the power system. Possible causes include incorrect handling of the instrument, use contrary to application, combination with an unsuitable power unit or saw template, interaction with other substances, insufficient identification marking, incorrect disposal, possible slipping during start-up, incorrect storage or transport, failure to eliminate a worn instrument and incorrect clamping in the power system.

Combination:
- Observe the instructions in the instructions for use enclosed with the power system.
- Set saw blades in motion inside the template before applying to the bone.
- The cutting width of the saw blades used must not exceed the thickness of the saw template. The saw blade must be able to move freely inside the saw template.

Preparation at the place of use/Storage and transport:
Store and transport the instruments at room temperature, in their original packaging away from dust, moisture and sunlight as well as sources of recontamination.

Storage:
Store the instruments at room temperature away from dust, moisture and sunlight as well as sources of recontamination.

Disposal:
Potentially infectious, pointed and sharp items have to be disposed of in non-pierceable and non-breakable disposable containers.

Procedure in case of serious adverse events:
All serious adverse events occurred in connection with the device have to be reported to the manufacturer and the competent national authority.

Liability:
Applies outside the USA: The user is responsible for checking the product prior to use to ensure that it is suitable for the intended purpose. The user is responsible for the use of the instruments. In case of contributory negligence by the user, Gebr. Bräsegger partially or totally declines liability for resulting damages, particularly if these arise from non-observance of our recommendations, misuse or warnings as well as inadvertent misuse by the user. Applies in the USA: ALL WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE EXPRESSLY DISCLAIMED. ALL IMPLIED OR CONSEQUENTIAL DAMAGES, INCLUDING, BUT NOT LIMITED TO, LOST PROFITS, INJURY TO OTHER PROPERTY, LOSS OF USE, LOSS OF OTHER COMMERCIAL LOSSES, ARE EXPRESSLY DISCLAIMED. IF, notwithstanding the foregoing, consequential and incidental damages cannot be excluded due to operation of law, such damages are expressly limited in amount to the purchase price of any defective products.

This medical device is marked as a disposable product. It is not intended for repeated use. Therefore, there are no recommendations to ensure safe reprocessing of the product.

Trademark Reference:
The third party trademarks used herein are trademarks of their respective owners.

(EN)	Symbols only ISO 15223-1: Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
(REF)	Catalogue number, ref. no. 5.1.6: Indicates the manufacturer's catalogue number so that the medical device can be identified
(LOT)	Batch code, ref. no. 5.1.5: Indicates the manufacturer's batch code so that the batch or lot can be identified
(STERILE)	Sterilized using irradiation, ref. no. 5.2.4: Indicates a medical device that has been sterilized using irradiation
(S)	Do not use if package is damaged, ref. no. 5.2.8: Indicates a medical device that should not be used if the package has been damaged or opened.
(M)	Manufacturer, ref. no. 5.1.1: Indicates the medical device manufacturer.
(D)	Date of manufacture, ref. no. 5.1.3: Indicates the date when the medical device was manufactured.
(D)	Distributor, ref. no. 5.1.9: Indicates the entity distributing the medical device into the locale
(U)	Use-by date, ref. no. 5.1.4: Indicates the date after which the medical device is not to be used.
(I)	Consult instructions for use, ref. no. 5.4.3: Indicates the need for the user to consult the instructions for use.
(S)	Single sterile barrier system, ref. no. 5.2.11: Indicates a single sterile barrier system.
(S)	Double sterile barrier system, ref. no. 5.2.12: Indicates two sterile barrier systems.
(S)	Single sterile barrier system with protective packaging inside, ref. no. 5.2.13: Indicates a single sterile barrier system with protective packaging inside.
(S)	Single sterile barrier system with protective packaging outside, ref. no. 5.2.14: Indicates a single sterile barrier system with protective packaging outside.
(S)	Do not re-use, ref. no. 5.4.2: Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.

Caution, ref. no. 5.4.4: Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.

Store in a dry place, ref. no. 5.3.4: Indicates that the medical device has to be protected from moisture.

Protect from sunlight, ref. no. 5.3.2: Indicates that the medical devices has to be protected from sources of light.

Contains hazardous substances, ref. no. 5.4.10: Indicates a medical device that contains substances that can be carcinogenic, mutagenic, reprotoxic (CMR), or substances with endocrine disrupting properties.

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